



## Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19

**Clinical Bulletin Number:** CB#010  
**Month/Year:** April 2020  
**Distribution:** Global  
**Date of implementation:** Immediately  
**Release version:** 1.0



This clinical bulletin provides clinicians treating COVID-19 patients with factors to consider when seeking methods to reduce the dispersion of aerosolized virus, and increase FiO<sub>2</sub> in the healthcare setting.

This low dispersion circuit is not the intended use for these devices, but emerging clinical practice in treating COVID-19 shows that it is a likely need given the ventilator shortage and the concern of infection risk to healthcare workers during the COVID-19 pandemic.

This circuit configuration applies to continuous positive airway pressure (CPAP), spontaneous (S), spontaneous timed (ST), pressure assist control (PAC), and timed (T) modes in Lumis, AirCurve, AirSense and S9 platforms. These devices are identified differently throughout global regions. **Section 6 describes the appropriate device name per region.**

Patient and device instructions for these devices should continue to be followed in addition to the information provided within this bulletin. **A summary of warnings and cautions are in section 5 of this bulletin.**

### 1. Circuit Setup

The low dispersion circuit is depicted below:




Figure 1: Low dispersion circuit setup

## Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19

The low dispersion circuit features, in succession: a Non Vented mask, an oxygen sideport connector, an Antibacterial/antiviral filter or Heat Moisture Exchange Filter, a ResMed Leak Valve and standard tubing.

Circuit component	Part description	Part number
Anti-bacterial (AB) / Anti-viral (AV) filter	Antibacterial filter – Adult, 50 pk	24966
Heat moisture exchanger filter (HMEF)	HME filter – Adult, 50 pk	24967
Oxygen port	Sideport oxygen connector	14915
Leak valve	ResMed leak valve ROW	24988



ResMed Leak Valve



14915 Oxygen Sideport Connector



26960 Oxygen Connector Port\*



Quattro FX NV



Quattro Air NV



AcuCare F1-0

**Figure 2. ResMed Circuit Ancillary Materials**

\*The 26960 oxygen connector port is also compatible with this circuit configuration



## Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19

### 2. Use of Low Dispersion Circuit and Supplemental Oxygen with CPAP and Bi-level Devices

ResMed has reviewed the risks associated with the use of this circuit configuration (with and without supplemental oxygen), and considers these risks to be clinically acceptable with the warnings and cautions listed in Section 5.

The assessment was performed on CPAP, S, ST, PAC and T modes in the following devices: S9 VPAP ST/ST-A, Lumis ST/ST-A, AirCurve 10 ST, S9 Elite (CPAP mode only), S9 AutoSet (CPAP mode only) and AirSense 10 Elite (CPAP mode only) and AirSense 10 AutoSet (CPAP mode only).

The Intelligent Volume-Assured Pressure Support (iVAPS) therapy mode was not assessed and is not recommended for this circuit configuration.

All mentioned devices do not comply to essential performance specifications required for a life support ventilator. The key differences between these devices and an invasive ventilator are:

- Lack of therapy alarms
- No internal battery
- No oxygen blender.

#### 2.1 CO<sub>2</sub> retention

Circuit accessories (eg, AB/AV filter, oxygen port, non-vented mask) placed between the leak port and the patient, contribute to increased dead-space which may lead to some rebreathing.

In order to minimize this, clinicians may consider switching to bi-level mode instead of CPAP, and/or increasing pressure support if patient CO<sub>2</sub> levels increase. See warnings and cautions in section 5.

#### 2.2 Mask Setting

When using this circuit configuration, it is recommended that the user select either the “nasal” or “full face” mask setting.

## Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19

### 3. Ventilator Performance Impact

Testing and clinical assessment by ResMed showed the ventilator performance was acceptable, but there may be a few impacts when Lumis, AirCurve, AirSense and S9 devices are used with a low dispersion circuit and with an oxygen flowrate of higher than 15L/min. **Also see warnings and cautions in section 5 of this bulletin.**



Device	Lumis ST and AirCurve ST	Lumis ST-A, AirCurve ST-A, S9 VPAP ST-A	AirSense 10 Elite and AirSense 10 AutoSet	S9 Elite, 29 AutoSet, S9 VPAP ST
Location	Hospital use only			
Oxygen	<ul style="list-style-type: none"> <li>Indicated for up to 15L/min</li> <li>May affect device performance with increasing O<sub>2</sub> levels</li> <li>Mitigate risk through external monitoring</li> </ul>			
Supervision	Monitoring with alarms strongly recommended			
Power	<ul style="list-style-type: none"> <li>Loss of power will stop ventilation</li> <li>No internal battery</li> </ul>			
Humidification or heated tubes	<ul style="list-style-type: none"> <li>Do not use integrated humidifier or heated tubes</li> <li>Heat-moisture exchanger (HME) may be used</li> </ul>			
Alarms	No integrated alarms, no alarm if ventilation stops. Notification on the screen will alert on high leak and blocked tube.	<ul style="list-style-type: none"> <li>Total power failure</li> <li>Blocked tube</li> <li>Disconnection</li> <li>Low pressure</li> <li>Low minute ventilation</li> </ul>	No integrated alarms, no alarm if ventilation stops. Notification on the screen will alert on high leak and blocked tube.	No integrated alarms, no alarm if ventilation stops. Notification on the screen will alert on high leak and blocked tube.
Dead space	AB filter, oxygen port, NV mask contribute to dead space, adjust pressure support to compensate.			
iVAPS	Do not use iVAPS mode as target volume accuracy will be degraded		iVAPS mode not available	iVAPS mode not available



## Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19


### 4. Supplemental oxygen

With increasing O<sub>2</sub> flow levels, triggering may be less sensitive to patient efforts and trigger delay may be increased. This may lead to asynchrony for very low patient flows, but can be improved by increasing the trigger sensitivity.

Pressure delivery and device reported values may be impacted by increasing O<sub>2</sub> flow levels, but not to the extent to impact patient safety or significantly reduce accuracy.

#### 4.1 FiO<sub>2</sub>

ResMed performed bench testing to provide users with guidance on the expected FiO<sub>2</sub> when using the recommended circuit. A compliant lung model (Resistance = 5 cmH<sub>2</sub>O, Compliance = 50 ml/cmH<sub>2</sub>O<sup>1</sup>) is used to achieve pressure waveforms and minute ventilations representative of COVID-19 patients. The breath rate is 20 bpm. CPAP: active patient model with 20 bpm. Bi-level: passive patient model with device backup rate set to 20 bpm.



<sup>1</sup> COVID-19 Does Not Lead to a “Typical” Acute Respiratory Distress Syndrome. Am J Respir Crit Care Med. 2020 Mar 30


ResMed Confidential Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way.



## Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19

**FiO<sub>2</sub> Concentration Output per O<sub>2</sub> Flow Rate, NV Full Face Mask, Unintentional Leak = 0 L/min**

Mode	CPAP	Bilevel	Bilevel	Bilevel
Pressure	10 cmH <sub>2</sub> O	15/5 cmH <sub>2</sub> O	16/8 cmH <sub>2</sub> O	20/10 cmH <sub>2</sub> O
Minute Ventilation	~7 L/min	~9 L/min	~7 L/min	~9 L/min
O <sub>2</sub> Flow Rate	FiO <sub>2</sub>			
5 L/min	54%	59%	45%	41%
10 L/min	78%	74%	68%	61%
15 L/min	94%	96%	94%	85%



**FiO<sub>2</sub> Concentration Output per O<sub>2</sub> Flow Rate, NV Full Face Mask, Unintentional Leak = 20 L/min**

Mode	CPAP	Bilevel	Bilevel	Bilevel
Pressure	10 cmH <sub>2</sub> O	15/5 cmH <sub>2</sub> O	16/8 cmH <sub>2</sub> O	20/10 cmH <sub>2</sub> O
Minute Ventilation	~7 L/min	~9 L/min	~7 L/min	~9 L/min
O <sub>2</sub> Flow Rate	FiO <sub>2</sub>			
5 L/min	47%	45%	45%	41%
10 L/min	70%	60%	66%	58%
15 L/min	85%	71%	79%	72%


ResMed Confidential Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way.



## Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19

When interpreting  $FiO_2$  values, it is important to keep in mind that various factors can influence the end results, both in bench testing and in the clinical application.

As can be seen in the graphs, both the **pressure settings and the resulting minute ventilation** create variation of at least 20%  $FiO_2$ . Additional variation may be caused by other confounding factors including but are not limited to: patient breath rate, inspiratory peak flows, unintentional<sup>2</sup> leak, and intentional leak of the ResMed leak valve.



Generally,  $FiO_2$  values will range between 50% and 80% at 10 L/min of  $O_2$  supply, and between 70% and 95% at 15 L/min of  $O_2$  supply, accounting for up to 20 L/min of unintentional mask leak.

<sup>2</sup> Unintentional leak refers to the air leak at the mask face interface, as opposed to intentional leak through the leak valve.



## Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19

### 4.2 Blocked tube alert

All devices contain a blocked tube alert/mitigation that may be falsely triggered with >15 L/min of supplemental O<sub>2</sub>. Devices will stop therapy if this mitigation occurs. If supplemental oxygen is added to the circuit at higher than 15 L/min, this might lead to false activation of 'blocked tube' alert, which leads to ventilation stop. Only Lumis ST-A, AirCurve ST-A and S9 VPAP ST-A will signal this alert with an audible alarm. The remaining listed devices do not contain audible alarms however, a notification on the screen will signal the blocked tube alert and lead to ventilation stop.

If using high flow levels of oxygen, external monitoring and supervision is recommended. If blocked tube alert occurs and ventilation stops, clinician may opt to increase pressure support, or decrease level of supplemental oxygen flow. Observe at least five minutes of the device operation during patient breath triggering when adding high flow rate of oxygen. See warnings and cautions on section 5.

### 4.3 Non-Vented mask alarm

An increased level of oxygen increases the likelihood of false triggering of the Non-Vented Mask alarm in Lumis ST-A, AirCurve ST-A and S9 VPAP ST-A. Clinicians should assess if opting to turn off the alarm in this circuit configuration would be appropriate and suit patient needs.





## Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19

### 5. Warnings and Cautions

When using the low dispersion circuit configuration, it is subject to the addition of the following:

**Warning:** This circuit, or the use of supplemental oxygen at the mask, may impact device performance and the effectiveness of alarms. When using this configuration with patients whose condition may deteriorate without adequate ventilation, use of independent patient monitoring equipment with alarms is strongly recommended.

**Warning:** Use of the device with this circuit configuration shall be limited to emergency use due to resource limitation in a hospital only.

**Warning:** The use of supplemental oxygen in excess of 15 L/min with a vented circuit with a ResMed Leak Valve may result in false triggering of the 'blocked tube' mitigation leading to ventilation stop<sup>3</sup>. The audible alarm is only present in Lumis ST-A, S9 VPAP ST-A and AirCurve 10 ST-A (US only) devices. Large supplemental oxygen levels should be verified as compatible by observing the patient for at least five minutes initially to check that a 'blocked tube' event is not triggered.

If the 'blocked tube' mitigation occurs and ventilation stops:


- Ensure the patient is cared for.
- Check the air tubing, remove any blockages and/or decrease the level of supplemental oxygen and/or increase the overall pressure setting.
- Press the dial to clear the message and then press Start/Stop to restart the device.
- Continue to monitor the patient.
- If the problem persists, stop using the device and contact your local ResMed dealer or ResMed office.

**Caution:** Circuit accessories (eg, AB/AV filter, oxygen port, non-vented mask) placed between the leak port and the patient contribute to increased dead-space which may lead to some rebreathing. Clinicians should consider minimizing that dead-space and adjusting device settings (eg, increasing pressure support/ventilation) to compensate for the increased dead-space if increase in patient CO<sub>2</sub> level is detected.

**Caution:** iVAPS mode is not recommended with this circuit configuration, especially when used with levels of entrained oxygen greater than 4 L/min. While the therapy continues to provide pressures within the clinician's set boundaries (PSmin and PSmax), the accuracy of control of volume (Target Va for iVAPS) might be degraded, potentially affecting the efficacy of the treatment.

---

<sup>3</sup> Higher levels of supplemental oxygen and resistance circuits result in lower average output flow being measured by the device. If this level is low enough for approximately two consecutive minutes, the device may falsely record a 'blocked tube' event.



**ResMed**

**Clinical Bulletin**

## **Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19**

**Note:** As currently stated in Lumis ST-A, AirCurve ST-A and S9 VPAP ST-A Clinical guides, the use of supplemental oxygen with a vented circuit may result in false triggering of the Non-Vented mask alarm – with increased level of oxygen up to 20 L/min, it increases the likelihood of false triggering of the Non-Vented mask alarm, therefore the clinician should assess if it is best to keep the alarm On or to turn it OFF.



## Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19

### 6. Device Regional Identification

The described circuit configuration applies to Lumis, AirCurve, AirSense and S9 platforms. The table below describes the appropriate device commercial name and function differentiated per region. Please note this does not apply to adaptive servo-ventilation (ASV), AutoSet (APAP), iVAPS modes.

Product Name	Regional Availability			Available Modes					Pressure Range* (cmH <sub>2</sub> O)
	AMR	EMEA	AGM	CPAP	S	T	ST	PAC	
S9 Elite		•	•	•					4-20
S9 AutoSet		•	•	•					4-20
S9 VPAP ST	•	•	•	•	•	•	•	•	ROW 2-25 AMR 3-25
S9 VPAP ST-A	•	•	•	•	•	•	•	•	ROW 2-30 AMR 3-30
AirSense 10 AutoSet	•	•	•	•					4-20
AirSense 10 Elite	•	•	•	•					4-20
AirCurve 10 ST	•	•		•	•	•	•		ROW 2-25 AMR 3-25
AirCurve 10 ST-A	•			•	•	•	•	•	3-30
Lumis 100/150 ST		•	•	•	•	•	•	•	2-25
Lumis 100/150 ST-A		•	•	•	•	•	•	•	2-30

Legend:


**AMR:** Americas

**EMEA:** Europe, Middle East, Africa

**AGM:** Asia Growth Markets

**ROW:** Rest of World

\*CPAP is always 4-20



**ResMed**

**Clinical Bulletin**

## **Supplemental Oxygen and Low Dispersion Circuit with CPAP and Bi-level Devices for COVID-19**

### **7. Multipatient use**

See the relevant device Clinical guide for Multipatient use details.

For further information, please contact your nearest ResMed office. See [ResMed.com](https://www.resmed.com) for details.

©2020 ResMed. C283653/01 2020-05



## Clinical bulletin

<b>Doc No.</b>	C283653
<b>Page</b>	13 of 13
<b>Rev</b>	1.0

<b>Rev</b>	<b>Change note</b>	<b>Document prepared by:</b>	<b>Document checked by:</b>
1.0	Obj. ref C283653	Bruno Sicre	Audrey Dauphin

## Clinical Bulletin CB#010

**Title of Clinical Bulletin: Supplemental Oxygen and Low Dispersion Circuit in CPAP and Bi-level devices for COVID-19.**

### 1. Details

<b>Brief Details:</b>	Clinical Bulletin describes the use of a low dispersion circuit with CPAP and bi-level devices (Lumis, AirCurve, AirSense and S9) in the context of addressing a shortage of dedicated ICU ventilators during the COVID-19 pandemic crisis.
<b>Style:</b>	If printed - (single-sided) on A4 80gsm white bond.
<b>Colour of printing:</b>	Colour
<b>Text:</b>	As shown.
<b>Manufacturer:</b>	External Contractor or In House Photocopy

### 2. Distribution

<b>By Medical Affairs to:</b>	(i) Regional product managers (ii) Product trainers (iii) Clinicians
-------------------------------	--